
KNEE EXTENSION

FITTING INSTRUCTIONS

It is best to apply the JAS Knee device sitting at the very edge of a seat. Once the device is applied choose any comfortable position (lying or sitting) to carry out your 30-minute JAS session. If you are using the JAS device in a seated position, remove shoe to allow foot to slide easily along the floor.

STEP 1: SET DRIVE ARMS IN CORRECT EXTENSION HOLE SETTINGS

IMPORTANT: Drive arms must be in same hole setting on both sides.


Fig. 1

Figure 1, 2, 3

- If necessary, place drive arms in Mid Range hole settings (Narrow or Wide to correspond with Tower setting).
- The Mid Range drive arm setting will be necessary for patients with significant ROM loss in knee extension.


Fig. 2




Fig. 3

STEP 2: PLACE DEVICE OVER LEG


- Align device over leg, centering tower over knee.

STEP 3: SECURE THIGH CUFF STRAPS

- Pull and secure inner thigh straps first, then outer thigh straps.

STEP 4: ADJUST DEVICE ANGLE TO MATCH END RANGE

- Extend knee to comfortable end range.
- Ratchet device so calf cuff makes solid contact with shin.

STEP 5: SECURE CALF CUFF

- Pull and secure calf straps comfortably snug.

NOTE: The calf cuff can be adjusted up or down for comfortable fit, by loosening the knob above cuff.

- Begin treatment protocol by rotating ratchet as outlined on back page.

RECOMMENDED

JAS[®] TREATMENT PROTOCOL

STEP 1: Rotate therapy knob until a gentle, pain-free stretch is felt. Level 2-3 on stretch intensity scale.

STEP 2: Hold stretch position for 5 minutes. Before turning knob, re-assess stretch level:
1. If stretch intensity has decreased, rotate knob until you feel a 2-3 level stretch again.
2. If stretch intensity has not changed, leave in same position.
3. If stretch intensity has increased, rotate knob in opposite direction until 2-3 level stretch is felt.

STEP 3: Hold stretch position for another 5 minutes.
Repeat steps 2 and 3 for a series of six 5-minute stretches (30-Minute Therapy Session).

IMPORTANT:

If you are initiating JAS therapy early post surgery, you may need to begin with shorter therapy sessions (one to three 5-minute stretches per session), due to post operative symptoms of swelling, pain and sensitivity to passive stretch. If so, gradually increase to 30-minute sessions over a two week period, as tolerated.

STEP 4: When 30-minute session is complete, turn knob in opposite direction until stretch is relieved, then remove device. You may feel some joint stiffness following your JAS session. Gently move your joint back and forth to “cool down” and relieve the stiffness.

IMPORTANT:

You must evaluate stretch sensation after each 5-minute stretch period, before adjusting the Range of Motion (ROM) on your JAS device. If stretch intensity has not decreased, it is incorrect to turn further into the ROM. You will not achieve a desired result of permanent tissue stretch.

Aggressive or painful levels of stretch will not be therapeutic, and will yield counterproductive results - increased stiffness, swelling or pain, and little or no gains in joint ROM.

NUMBER OF SESSIONS PER DAY

Three 30-minute sessions, per direction, per day, are suggested for optimal results. Begin with one 30-minute session per day and add sessions every 3-5 days as tolerated up to 3 sessions per direction per day.

WARNING: JAS SPS devices are NOT designed to be worn overnight or for prolonged, continuous periods of time. Joint Active Systems does not recommend or endorse a therapy protocol that includes extended wear times of 4 + hour sessions as are commonly recommended with dynamic splints or turnbuckle splints. This is due to the known increased risk of skin irritation and / or breakdown associated with any prolonged splinting regime.

Wait 45-60 minutes between each JAS treatment session.

If you experience a change in swelling or irritation to your underlying condition when using the JAS device, discontinue use and contact your physician.

TECHNICAL SUPPORT

Call JAS toll free at **(800) 879-0117** for technical assistance and any questions regarding your JAS device. **This is a rental device.** At the end of your treatment it is **your responsibility** to contact JAS to terminate billing and arrange for device return.

JAS Representative

Contact Number

Joint Active Systems, Inc.

2600 South Raney • Effingham, IL 62401

TEL: (217) 342-3412 or (800) 879-0117

Email: info@jointactivesystems.com


www.jointactivesystems.com

Covered by one or more US patents. Other patents pending.

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